

# Next: QuickMIC® GN02

## Expanded Gram-negative AST

### for Highly Resistant Pathogens

Rapid MIC in 2-4 Hours Directly From Positive Blood Cultures

The data presented are based on preliminary findings<sup>1</sup>. Further studies and additional data are required to fully validate the performance of the final QuickMIC® GN02 panel configuration.

### QuickMIC® GN02 in development

**QuickMIC GN02** will broaden Gram-negative antibiotic coverage with advanced, resistance-targeted agents designed to overcome key  $\beta$ -lactam resistance mechanisms. This will enable rapid phenotypic AST directly from positive blood cultures, supporting timely and effective treatment of infections caused by multidrug-resistant Gram-negative pathogens.

### Expanded antibiotic coverage

Designed for today's antimicrobial resistance challenges, QuickMIC GN02 is focused on newer antibiotics and inhibitor combinations, including cefiderocol, ceftolozane/tazobactam, meropenem/vaborbactam, and aztreonam/avibactam.

### Cefiderocol spotlight

**Cefiderocol** is a siderophore cephalosporin used as a last-resort therapy against resistant Gram-negative pathogens. Preliminary QuickMIC data suggests **rapid phenotypic AST** may be feasible directly from positive blood culture, with cefiderocol MIC results available in ~3.5 hours.

#### From testing challenge to rapid results


- ✓ Testing cefiderocol using conventional automated AST systems is challenging because the antibiotic requires non-standard testing conditions with specialised iron-depleted media.
- ✓ Preliminary data indicates that QuickMIC allows testing without the need for specialised iron-depleted media, supporting susceptibility testing for this critical therapy.



#### QuickMIC® GN02:

- ✓ Amikacin
- ✓ Aztreonam/avibactam
- ✓ Cefepime
- ✓ Cefiderocol
- ✓ Ceftazidime
- ✓ Ceftazidime/avibactam
- ✓ Ceftolozane/tazobactam
- ✓ Ciprofloxacin
- ✓ Colistin
- ✓ Meropenem
- ✓ Meropenem/vaborbactam
- ✓ Piperacillin/tazobactam

#### Cefiderocol results from performance testing



QuickMIC® data shows >90% Categorical Agreement for cefiderocol on a resistant challenge dataset exceeding 20 carbapenemase-producing E.coli<sup>1</sup>.

<sup>1</sup>Åman A, Chatzopoulou M, Öhrn H, Faxén L, Johansson C, Olsson A, Lagerbäck P, Tängdén T, Malmberg C. Rapid susceptibility testing of cefiderocol directly from positive blood cultures using a microfluidic agar-diffusion. ESCMID Global 2025, P1550, April 2025.

QuickMIC® GN02 supports susceptibility testing of newer antibiotics and inhibitor combinations for Gram-negative pathogens in clinical microbiology laboratories.